


pH Metry




Contents

- Aim
- Introduction
- Working Principle of instrument
- Materials and Method for sample analysis
- Applications



Aim: To study about pH Meter.



Introduction

- The pH meter was invented in 1934 by the American chemist Arnold O. Beckman (1900-2004) to measure the sourness of lemons.
- A simple and speedy device to measure the acidity and alkalinity of a fluid

➤ pH (Potential of Hydrogen):

The logarithm of the reciprocal of hydrogen-ion concentration in gram atoms per litre; provides a measure on a scale from 0 to 14 of the acidity or alkalinity of a solution (where 7 is neutral and greater than 7 is more basic and less than 7 is more acidic)

The formal definition of pH is the negative logarithm of the hydrogen ion activity.

$$\text{pH} = -\log[\text{H}^+]$$

Working Principle of instrument

- When a pair of electrodes or a combined electrode(glass & calomel electrode) is dipped in an aqueous solution, a potential is developed across the thin glass of the bulb. The e.m.f of complete cell(E) formed by linking of this two electrodes at a given solution temperature is therefore,


$$E = E_{\text{ref}} - E_{\text{glass}}$$

E_{ref} = The potential of stable calomel electrode which at normal room temp is +0.250V.

E_{glass} = The potential of glass electrode which depends on the pH of the solution under test.

A pH meter acts as a volt meter that measures the electrical potential difference between a pH electrode and a reference electrode and displays the result in terms of the pH value of the solution in which they are immersed.

An electrical potential develops when one liquid is brought into contact with another one ,but a membrane is needed to keep such liquid apart.




Materials and Method for sample analysis

➤ **A pH measurement system consists of**

- i) 2 probes
- ii) Buffer powder/capsules for calibration
- iii) High input meter
- iv) Distilled water

- The pH measuring electrode is a hydrogen ion sensitive glass bulb.
- The reference electrode output does not vary with the activity of the hydrogen ion.

Types of Electrode:




a. Calomel electrode: It consists of a glass tube containing saturated KCl connected to platinum wires through mercury (mercury chloride)

b. Glass electrode: It consists of a thin bulb of special glass blown at the end of the glass tube, and the bulb is filled with dilute acid for example decinormal HCl acid connected to a silver chloride electrode

c. Reference junction: typically made of porous ceramic or porous teflon



Combined Electrode:

- It is easier to handle one electrode instead of two.
- The indicating glass electrode and the reference electrode are simply built into a single physical entity.
- This helps to ensure that the two electrodes have the same temperature during operation.




Temperature and Buffers

- Temperature compensation is contained within the instrument because pH electrodes are temperature sensitive.
- Temperature compensation only corrects for the change in the output of the electrode, not for the change in the actual solution.
- Buffers are solutions that have constant pH values and the ability to resist changes in pH.
- They are used to calibrate the pH meter.




Method of Using

1. The 2 probes should connect to pH reading meter
2. Connect the power pack of 230v to the pH meter
3. Dip the electrode in a standard solution of pH and set the temperature and take reading [The 1st probe should be calibrate with different buffer solutions for rule out the errors of the electrode (buffer solution : buffer powder + distil water) Then should check the solution temperature with 2nd probe]
4. Remove buffers wash and wipe the electrode and dip it in the desirable solution and take the reading.




Applications:

1. Hospitals
2. Universities
3. Pharmaceutical companies
4. Research institutes or laboratories
5. Measure soil pH for healthier plants.




LOUD SPEAKER

- When things shake about, or vibrate, they make the [sounds](#) we can hear in the world around us.
- At the front of a loudspeaker, there is a fabric, [plastic](#), [paper](#), or lightweight metal **cone** (sometimes called a **diaphragm**) not unlike a drum skin (colored gray in our picture).
- the inner part is fixed to an [iron](#) coil (sometimes called the **voice coil**, colored orange in the diagram) that sits just in front of a permanent [magnet](#) (sometimes called the **field magnet**, and colored yellow).
- When you hook up the loudspeaker to a stereo, electrical signals feed through the speaker cables (red) into the coil.
- This turns the coil into a temporary magnet or **electromagnet**. As the [electricity](#) flows back and forth in the cables, the electromagnet either attracts or repels the permanent magnet.
- This moves the coil back and forward, pulling and pushing the loudspeaker cone. Like a drum skin vibrating back and forth, the moving cone pumps sounds out into the air.



Microphones

- When you speak, **sound waves** created by your voice carry energy toward the microphone. Remember that sound we can hear is energy carried by vibrations in the air.
- Inside the microphone, the **diaphragm** (much smaller than you'd find in a loudspeaker and usually made of very thin [plastic](#)) moves back and forth when the sound waves hit it.
- The **coil**, attached to the diaphragm, moves back and forth as well.
- The **permanent magnet** produces a [magnetic field](#) that cuts through the coil. As the coil moves back and forth through the magnetic field, an [electric current](#) flows through it.
- The **electric current** flows out from the microphone to an amplifier or sound recording device. Hey presto, you've converted your original [sound](#) into electricity! By using this current to drive sound recording equipment, you can effectively store the sound forever more. Or you could [amplify](#) (boost the size of) the current and then feed it into a loudspeaker, turning the electricity back into much louder sound. That's how PA (personal address) systems, [electric guitar](#) amplifiers, and rock concert amplifiers work.



Thank you...

